

MASTERPULS® MP50 »ultra«



- Weight main device: 9.5 kg
- Pressure: 1.0 – 4.0 bar
- Frequency: 1 – 15 Hz
- Vibration therapy, V-ACTOR®: 21 Hz
- Built-in standard air compressor

MASTERPULS® MP100 »ultra«



- Weight main device: 10.5 kg
- Pressure: 1.0 – 5.0 bar
- Frequency: 1 – 21 Hz
- Vibration therapy, V-ACTOR®: 31 Hz
- Extremely silent »Air Power« drive

General data MP50/100



- Housing dimensions (W x H x D):
426 x 144 x 340 mm
- Weight R-SW handpiece: 480 g
- Weight V-ACTOR® handpiece: 370 g
- Two handpiece connectors
- Class IIa devices
- Online access to the »ICE Shock Wave Portal«
- Transmitter case (optional)
- Trolley (optional)

»ICE Shock Wave Portal«



- Selection of indications
- Treatment instructions
- Treatment videos

